For the Quantitative Measurement of free T₄ in Human Serum

The concentration of free T₄ (FT₄) in a given specimen determined with assays from different manufacturers can vary due to differences in assay methods and reagent specificity. The results reported by the laboratory to the physician must include the identity of the free T₄ assay method used. Values obtained with different assay methods should not be used interchangeably.

CAUTION: United States Federal law restricts this device to sale and distribution by or on the order of a physician, or to a clinical laboratory; and use is restricted to, by or on the order of a physician.

INTENDED USE

The FastPack® IP and FastPack® free T₄ Immunoassay is a paramagnetic particle immunoassay for the *in vitro* quantitative determination of FT₄ in human serum. The FastPack® IP and FastPack® free T₄ Immunoassay is designed for use with the FastPack® IP System.

SUMMARY

The thyroid hormone, thyroxine (T₄) is secreted from the thyroid gland located in the neck. Thyroid hormone is indispensable for growth, development, and sexual maturation. T₄ effects metabolic activity in many tissues leading to increased consumption of oxygen, and stimulation of mitochondrial respiration¹. Additionally, T₄ action includes stimulation of heart rate and heart contraction, stimulation of protein synthesis and carbohydrate metabolism, increase in the synthesis and degradation of cholesterol and triglycerides, increase in vitamin requirements, and enhancement of sensitivity of β-adrenergic receptors to catecholamines. These effects are usually magnified in patients with hyperthyroidism and minimized in patients with hypothyroidism².

T₄ in the circulation is reversibly and almost completely bound to carrier proteins (T₄-binding globulin, T₄-binding prealbumin, and albumin). Only about 0.03 % of the T₄ is unbound or free. All of the biological activity is associated with the free form of the hormone (free T₄). The large concentration of protein-bound hormone acts as a reservoir of T₄, maintaining a constant concentration of free T₄. Large variations in the concentration of these proteins, and total T₄ levels even in normal (euthyroid) individuals is temporized by having the active, free-form of the hormone as a small fraction of the total. Hence measurement of the free T₄ level is important in the diagnosis of hormone excess or deficient states.

Measuring the free T₄ has traditionally been accomplished by indirect means. Using an assay for the total T₄ coupled with a Thyroxine Uptake (TU), a Free Thyroxine Index (FTI)³, which is indirectly proportional to the FT₄ value, can be obtained. Direct determinations use equilibrium dialysis, ultrafiltration, RIA, and solid-phase EIAs are used to determine the FT₄ levels.

The FastPack® IP and FastPack® free T₄ Immunoassay is based on a competitive immunoassay principle: free T₄ in the sample (control or patient) competes with T₃ covalently coupled to a solid phase paramagnetic particle for an alkaline phosphatase-labeled anti-T₄ monoclonal antibody. The amount of labeled-anti-T₄ antibody bound to the magnetic particle is inversely proportional to the concentration of FT₄ in the sample.

TEST PRINCIPLE

The FastPack® IP and FastPack® free T₄ Immunoassay is a competitive chemiluminescence assays.

- Primary incubation: Antibody solution (a buffer solution containing a monoclonal anti-thyroxine (T₄) antibody labeled with alkaline phosphatase)[100 μL] reacts with Free T₄ from the patient's sample, control, or calibrator [100 μL].
- Secondary incubation: The reaction mixture is added to streptavidin coated paramagnetic particles with tri-iodo thyronine (T₃)-biotin on the particle surface. During this incubation, the T₃-coated beads compete with the sample free T₄.
- Removal of unbound materials: The paramagnetic particles are washed with wash buffer [0.2 mL/wash] to remove unbound materials.
- Substrate addition and detection: Chemiluminogenic substrate [140 μL] is added to the solid-phase bound complex and results in "glow" chemiluminescence, which is measured using the FastPack® IP Analyzer.
- The amount of bound labeled-antibody is inversely proportional to the concentration of FT₄ in the sample.

REAGENTS – Content and Concentration

Each FastPack® IP carton contains:

- 30 FastPack® IPs

Each FastPack® carton contains:

- 50 FastPacks

Each FastPack® IP or FastPack® Contains:

- Paramagnetic Particles, 150 µL
- Streptavidin-coated paramagnetic particles with bound T₃-biotin in buffer containing 0.1% sodium azide as a preservative.
- T₄ Antibody Solution, 100 µL
- Antibody solution containing mouse monoclonal anti-T₄ antibody labeled with alkaline phosphatase in a protein matrix containing Proclin® 150 as preservative.
- Wash Buffer, 2.0 mL
- Tris buffer containing surfactants.
- Substrate, 140 µL
- ImmGlow™: Indoxyl-3-phosphate and lucigenin in buffer containing preservatives.

Materials required but not provided

- FastPack® IP System
- FastPack® free T₄ Calibrator Kit – Cat. No. 25000016
- FastPack® free T₄ Control kit – Cat No. 25000027

WARNINGS AND PRECAUTIONS

- For *in vitro* diagnostic use only.
- Do not pipette by mouth.
- Do not eat, drink or smoke in designated work areas.
- Wash hands thoroughly after handling specimen.
- HAMA Interference: some individuals have antibodies to mouse protein (HAMA), which can cause interference in immunoassays that employ antibodies derived from mice⁶.
- FastPack® IP and FastPack® reagents are stable until the expiration date on the label when stored and handled as directed. Do not use FastPack® IP and FastPack® reagents beyond the expiration date.
- Discard used FastPack® IPs into a Biohazard container.
- The components containing Sodium Azide are classified per applicable European Economic Community (EEC) Directive as: Very toxic and dangerous to the environment (T+ N). The following are the appropriate Risk (R) and Safety (S) phrases:

R28	Very toxic if swallowed
R32	Contact with acids liberate very toxic gas
R50/53	Very toxic to aquatic organisms may cause long-term adverse effects in the aquatic environment.
S28	After contact with skin, wash immediately with plenty of soap suds.
S45	In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible).
S60	This material and its container must be disposed of a hazardous waste.
S61	Avoid release to the environment. Refer to special instructions/safety data sheets.

- The components containing Proclin 150 are classified per applicable European Economic Community (EEC) Directives as: Corrosive (C). the following are appropriate Risk (R) and Safety (S) phrases for Proclin 150:

R34	Causes burns.
R43	May cause sensitization by skin contact.
S26	In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
S36/37/39	Wear suitable protective clothing, gloves, and eye/face protection.
S45	In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible)

STORAGE INSTRUCTIONS

Store at 2 – 8 °C. Protect from light.

SPECIMEN COLLECTION/PREPARATION

1. Serum samples can be used for the FastPack® IP and FastPack® free T₄ Immunoassay.
2. The National Committee for Clinical Laboratory Standards (NCCLS) provides recommendations for handling, processing, and storing blood.^{4,5}
3. Collect all blood samples observing routine precautions for venipuncture.
4. For serum samples:
 - Collect samples in red top tubes with no gel separator.
 - Ensure that complete clot formation has occurred before centrifugation. Some specimens may exhibit increased clotting time, especially those from patients receiving anticoagulant or thrombolytic therapy.
 - Serum should be separated from the clot and centrifuged within 3 hours from time of collection.

- Remove serum from the cells prior to storage at 2–8° C.
 - If not tested within 24 hours, the sample should be frozen at –20° C or colder.
5. Do not freeze samples (–20 °C) for more than two months.
 6. Frozen samples should be thawed completely and mixed by gentle inversion prior to use.
 7. Samples should be free of Fibrin, Red Blood Cells, or other particulate material for optimal results.
 8. Serum samples showing turbidity and particulate matter should be centrifuged prior to use.
 9. Ensure the samples are free of bubbles.
 10. Human samples should be handled in accordance with the OSHA standard on Blood borne Pathogens.⁷

ASSAY PROCEDURE

See the FastPack® IP System Procedure Manual for information on operating the FastPack® IP System.

INSTRUMENTATION

FastPack® IP System

DETAILS OF CALIBRATION

During the FastPack® IP and FastPack® production process, Qualigen generates a master standard curve and places this information in the barcode of each FastPack® IP and FastPack® label, where it can be read by the FastPack® IP Analyzer during the testing sequence. The FastPack® IP Analyzer must be calibrated by the user to ensure that it is properly adjusted for the particular lot of FastPacks that is being used. Separate calibrations must be run for each type of test, i.e. Total PSA, Testosterone, or free T₄. The frequency of calibration varies for each test type. For the FastPack® IP and FastPack® free T₄ Immunoassay, the FastPack® IP Analyzer must be calibrated once every 14 days or whenever a new lot of free T₄ FastPacks are to be used.

Whenever the user performs an initial calibration for a particular lot of FastPacks or uses a new lot of calibrator, 2 FastPacks must be run for calibration (duplicates). Whenever recalibration is performed with the same lot of FastPacks and calibrator, 2 FastPacks must be run for calibration. See FastPack® IP System Procedure Manual for “Running a Calibration”.

Use FastPack® free T₄ Calibrator Kit – Cat. No. 25000016

RESULTS

The FastPack® IP Analyzer uses the information from the barcode to construct a lookup table of x,y values that represent the standard curve and estimates the concentration of unknown samples by linear interpolation.

QUALITY CONTROL

Quality control materials simulate real specimens and are essential for monitoring the system performance of assays. Good Laboratory Practices (GLP) includes the use of control specimens to ensure that all reagents and protocols are performing properly. Commercial control products (such as Bio-Rad Liqui-chek® Immunoassay Plus Control, level 2 and 3) are available which contain FT₄. See FastPack® IP System Procedure Manual for “Running Controls”.

LIMITATION OF PROCEDURE

- Specimens can be measured within the reportable range of the analytical sensitivity and the upper end of the calibration range, 6 ng/dL.
- Samples >6 ng/dL should be reported as such or re-run using another method. Dilution of out of ranges results is not recommended.
- Specimen from patients who have received preparations of mouse monoclonal antibodies for diagnosis or therapy may contain human anti-mouse antibodies (HAMA). Such specimens may show either falsely elevated or depressed values when tested with assay kits employing mouse monoclonal antibodies.
- Heterophilic antibodies in a sample have the potential to cause interference in immunoassay systems. Infrequently, FT₄ levels may appear elevated due to heterophilic antibodies present in the patient’s serum or to nonspecific protein binding. If the FT₄ level is inconsistent with clinical evidence, additional FT₄ testing is suggested to confirm the result.
- For diagnostic purposes, the FastPack® IP and FastPack® free T₄ Immunoassay should always be assessed in conjunction with the patient’s medical history, clinical examination and other findings.

EXPECTED RANGE

Serum samples were obtained from 207 individuals (99 male and 108 female). Samples were obtained from normal, healthy blood donors without any clinically abnormal indications. FT₄ levels were determined using the FastPack® IP and FastPack® free T₄ Immunoassay in conjunction with the FastPack® IP Analyzer in order to establish the FT₄ concentration in the normal population. From this data the expected normal range for the FastPack® IP and FastPack® free T₄ Immunoassay is 0.5 – 1.4 ng/dL. The expected range reflects the 95% Confidence Interval of the donor population of this study group. Each laboratory should determine their own reference range appropriate for their population.

SPECIFIC PERFORMANCE CHARACTERISTICS

Precision

The reproducibility of the FT₄ assay was measured by 2 patient samples in duplicate (n=54, for each sample) over ten days using three analyzers and two lots of reagents. Precision is expressed as the Coefficient of Variation.

Sample	Mean (ng/dL)	Between Run CV%	Between Analyzer CV%	Between Reagent Lot CV%	TOTAL CV%
1	1.43	7.4	2.8	11.1	11.5
2	3.51	6.2	1.8	4.8	7.2

Method Comparison

Clinical samples were used to compare the serum values obtained using the FastPack[®] IP and FastPack[®] FT₄ method and the serum value obtained using the Abbott IMx Free T₄ method. The values were evaluated for agreement using Deming regression analysis. Pearson correlation was used to test association between the values.

n	Range of Observation (ng/dL)	Intercept (ng/dL)	Slope	r
131	0.41 to 5.93	-0.169	1.07	0.95

Interfering Substances

Interfering substances were added to serum samples containing known amounts of FT₄. The value obtained for the serum with each interfering substance was compared to the value obtained for the serum without the interfering substance. These compounds did not show interference at the levels indicated.

Test Compound	Test Concentration
Bilirubin	40 mg/dL
Hemoglobin	1,000 mg/dL
Triglycerides	1,000 mg/dL

Analytical Sensitivity

The analytical sensitivity of FT₄ is 0.4 ng FT₄/dL. This value was interpolated from the dose-response curve that was minus two standard deviations below the average of the zero calibrator signal (n >10/analyzer/lot).

Specificity

Cross-reactivity of the antibody used in the FastPack[®] IP and FastPack[®] Free T₄ assay was shown to be < 1% for L-T₃.

REFERENCES

- ¹ Braverman, LE, Utigen, RD, Eds: Werner and Ingbar's The Thyroid. A Fundamental and Clinical Text. 7th ed., Philadelphia, Lippincott-Raven, 1996.
- ² Burtis, CA, Ashwood, ER, Eds: Tietz Textbook of Clinical Chemistry, 3rd ed., New York, WB Saunders Company, p. 1492, 1999.
- ³ Horworth, PJN, Ward RL, J. Clin.Pathol., 25, 259, 1972.
- ⁴ Approved Standard – procedures for the collection of diagnostic blood specimens by venipuncture. 5th Edition: H3-A5: 23(32), 2003, National Committee for Clinical Laboratory Standards (NCCLS)
- ⁵ Approved guideline – procedures for the handling and processing of blood specimens, H18-A2;19(21), 1999. National Committee for Clinical Laboratory Standards (NCCLS).
- ⁶ Schroff, RJ, Foon, KA, et.al.: Human anti-mouse immunoglobulin responses in patients receiving monoclonal antibody therapy. Cancer Res, 45:879 – 885, 1985
- ⁷ US Department of Labor, Occupational Safety and Health Administration, 29CFR Part 1910.1030, Occupational Exposure to Bloodborne Pathogens; Final rule. Federal Register 1991; 56(235): 64175-82.

Qualigen, Inc.
Carlsbad, CA 92011 USA
Technical Support
(760) 918-9165
(877) 709-2169

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© 2000 Qualigen, Inc. All rights reserved. Qualigen and FastPack are trademarks or registered trademarks of Qualigen, Inc.
All other trademarks are the property of their respective owners.

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