For use in conjunction with the FastPack® IP and FastPack® TSH Immunoassay and FastPack® IP System Analyzer

**CAUTION: United States Federal law restricts this device to sale and distribution by or on the order of a physician, or to a clinical laboratory; and use is restricted to, by or on the order of a physician.**

### INTENDED USE

The FastPack® TSH Controls are assayed quality control materials for the verification of the accuracy and precision of the FastPack® IP System when used for the quantitative determination of TSH in human plasma.

### SUMMARY AND PRINCIPLE

The use of control material is indicated as an objective assessment of the precision of methods and techniques in use. Two levels of control are available to allow performance monitoring within the clinical range.

### PRODUCT INFORMATION

- Provided ready to use.
- Mix contents by gently inverting before use. Avoid bubble formation.
- Controls: 5 mL/vial. Liquid. Contains components of human origin prepared in a Tris buffer solution with protein stabilizers to yield predetermined concentrations.
  - Control 1
  - Control 2
- Preservative: 0.1% sodium azide

### WARNINGS AND PRECAUTIONS

- **For *In Vitro* diagnostic use only.**
- Do not pipette by mouth.
- Do not eat, drink, or smoke in designated work areas.
- Do not mix controls from different lots.
- After opening, controls are stable until the expiration date on the label when stored and handled as directed. Do not use controls beyond the expiration date. Control stability is based on a closed vial stability study.
- Avoid microbial contamination of reagent when removing aliquots from the bottles.
- Refer to the FastPack® IP Procedure Manual for calibration procedures.
- Discard unused or expired control material, in stoppered vial, into a Biohazard container.
- Reagents in this kit contain sodium azide as a preservative, which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of such reagents, always flush with large volumes of water to prevent azide build up.
- The components containing sodium azide are classified per applicable European Economic Community (EEC) Directives as: Very toxic and dangerous to the environment (T+ N). The following are appropriate Risk (R) and Safety (S) phrases for sodium azide:
  - R28 Very toxic if swallowed.
  - R32 Contact with acids liberates very toxic gas.
  - R50/53 Very toxic to aquatic organisms. May cause long-term adverse effects in the aquatic environment.
  - S28 After contact with skin, wash immediately with plenty of soap-suds.
  - S45 In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible).
  - S60 This material and its container must be disposed of as hazardous waste.
  - S61 Avoid release to the environment. Refer to special instructions/safety data sheets.
- **Human source material. The antigens used in the preparation are potentially infectious and should be handled according to universal precautions and good clinical laboratory practices. Where appropriate, the donors were screened for HIV, HBV and HCV using FDA approved tests.**

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Qualigen, Inc.