For use in conjunction with the FastPack® IP and FastPack® TSH Immunoassay and FastPack® IP System Analyzer

CAUTION: United States Federal law restricts this device to sale and distribution by or on the order of a physician, or to a clinical laboratory; and use is restricted to, by or on the order of a physician.

INTENDED USE

The FastPack® TSH Calibrator is intended to calibrate the FastPack® IP System when used for the quantitative determination of TSH in human serum and plasma.

SUMMARY AND PRINCIPLE

Quantitative assay calibration is the process by which a set of samples with known analyte concentrations (i.e. assay calibrators) is tested to measure the response. The mathematical relationship between the measured responses and the known analyte concentrations establishes the calibration curve. The FastPack® IP System analyzer must be calibrated by the user to ensure that it is properly adjusted for the particular lot of FastPacks that is being used. The FastPack® TSH Calibrator is used for this purpose.

PRODUCT INFORMATION

- Provided ready to use.
- Mix contents by gently inverting before use. Avoid bubble formation.
- Calibrator: 5 mL/vial. Liquid. Contains components of human origin prepared in a Tris buffer solution with protein stabilizers to yield predetermined concentrations.
Calibrator B
- Calibration Card: 1
- Preservative: 0.1% sodium azide

WARNINGS AND PRECAUTIONS

- **For *In Vitro* diagnostic use only.**
- Do not pipette by mouth.
- Do not eat, drink, or smoke in designated work areas.
- Do not mix calibrators from different lots.
- After opening, calibrators are stable until the expiration date on the label when stored and handled as directed. Do not use calibrators beyond the expiration date. Calibrator stability is based on a closed vial stability study.
- Avoid microbial contamination of reagent when removing aliquots from the bottles.
- Refer to the FastPack® IP Procedure Manual for calibration procedures.
- Discard unused or expired calibrator material, in stoppered vial, into a Biohazard container.
- Reagents in this kit contain sodium azide as a preservative, which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of such reagents, always flush with large volumes of water to prevent azide build up.
- The components containing sodium azide are classified per applicable European Economic Community (EEC) Directives as: Very toxic and dangerous to the environment (T+ N). The following are appropriate Risk (R) and Safety (S) phrases for sodium azide:
 - R28 Very toxic if swallowed.
 - R32 Contact with acids liberates very toxic gas.
 - R50/53 Very toxic to aquatic organisms. May cause long-term adverse effects in the aquatic environment.
 - S28 After contact with skin, wash immediately with plenty of soap-suds.
 - S45 In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible).
 - S60 This material and its container must be disposed of as hazardous waste.
 - S61 Avoid release to the environment. Refer to special instructions/safety data sheets.
- **Human source material. The antigens used in the preparation are potentially infectious and should be handled according to universal precautions and good clinical laboratory practices. Where appropriate, the donors were screened for HIV, HBV and HCV using FDA approved tests.**

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© 2000 Qualigen, Inc. All rights reserved. Qualigen and FastPack are trademarks or registered trademarks of Qualigen, Inc.
All other trademarks are the property of their respective owners.

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