For use in conjunction with the FastPack® IP and FastPack® PSA Immunoassay and the FastPack® IP System Analyzer

**CAUTION: United States Federal law restricts this device to sale and distribution by or on the order of a physician, or to a clinical laboratory; and use is restricted to, by or on the order of a physician.**

### INTENDED USE

The FastPack® PSA Controls are assayed quality control materials for the verification of the accuracy and precision of the FastPack® IP System when used for the quantitative determination of PSA in human serum and plasma.

### SUMMARY AND PRINCIPLE

The use of control material is indicated as an objective assessment of the precision of methods and techniques in use. Two levels of control are available to allow performance monitoring within the clinical range.

### PRODUCT INFORMATION

FastPack® PSA Controls - Cat. No. 25000003

Note: The control range card for concentrations of the controls are located in the box.

- Provided ready to use.
- Mix contents by gently inverting before use. Avoid bubble formation.
- Controls: 5 mL/vial. Liquid. Contain components of human origin prepared in a Tris Buffer with protein stabilizers to yield predetermined concentrations.
  - Control 1
  - Control 2
- Preservatives: 0.1% Sodium Azide

### WARNINGS AND PRECAUTIONS

- **For *in vitro* diagnostic use only.**
- Do not pipette by mouth.
- Do not eat, drink, or smoke in designated work areas.
- Do not mix controls from different lots.
- After opening, controls are stable until the expiration date on the label when stored and handled as directed. Do not use controls beyond the expiration date.
- Avoid microbial contamination of reagent when removing aliquots from the bottles.
- Refer to the FastPack® IP Procedure Manual for additional information.
- Discard unused or expired control material, in stoppered vial, into a Biohazard container.
- Reagents in this kit contain sodium azide as a preservative, which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of such reagents, always flush with large volumes of water to prevent azide build up.
- The components containing Sodium Azide are classified per applicable European Economic Community (EEC) Directives as: Very toxic and dangerous to the environment (T+ N). The following are the appropriate Risk (R) and Safety (S) phrases:
  - R28 Very toxic if swallowed.
  - R32 Contact with acids liberates very toxic gas.
  - R50/53 Very toxic to aquatic organisms may cause long-term adverse effects in the aquatic environment.
  - S28 After contact with skin, wash immediately with plenty of soap-suds.
  - S45 In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible).
  - S60 This material and its container must be disposed of as hazardous waste.
  - S61 Avoid release to the environment. Refer to special instructions/safety data sheets.
- **Human source material. Treat as potentially infectious.**

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© 2000 Qualigen, Inc. All rights reserved. Qualigen and FastPack are trademarks or registered trademarks of Qualigen, Inc.  
All other trademarks are the property of their respective owners.

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