

The FastPack[®] IP System

PRODUCT SPECIFICATIONS

TECHNOLOGY:
Paramagnetic Particle
Chemiluminescence Immunoassay

OPERATING TEMPERATURE:
15°C (59°F) to 32°C (90°F)

OPERATING HUMIDITY:
10% to 80% Relative Humidity

DISPLAY OUTPUT:
2-line Alphanumeric Display

DATA OUTPUT:
RS 232, 9600 Baud, ASCII Text


FIELD SERVICE OPTIONS:
ROM Card Firmware Upgrade

ELECTRICAL REQUIREMENTS:
100-240 ± 10% VAC, 1.9A

POWER CONSUMPTION:
100 Watts (Maximum)

SIZE:
13"W x 9"H x 12"D (33 x 23 x 30 cm)

WEIGHT:
27 lbs (12.25kg)



Qualigen[®]

www.qualigeninc.com